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**MANUFACTURER IDENTIFICATION:**

Type: FFP1 NR; FFP2 NR; FFP3 NR

NR is to be worn for a single shift only

D: Meet the Dolomite clogging test

CE Standard considered : EN149:2001+A1:2009

CE Marking: CE 2797

#### Technical information

Tests	Type	FFP1 NR	FFP2 NR	FFP3 NR
Total Inward Leakage		<22%	<8%	<2%
Penetration	S	<20%	S	<6%
	L	<20%	L	<6%
Carbon dioxide content of inhalation air (%)		<1%	<1%	<1%
Inhalation resistance 30 l/min (mbar)		<0.6 mbar	<0.7 mbar	<1 mbar
Inhalation resistance 95 l/min (mbar)		<2.1 mbar	<2.4 mbar	<3 mbar
Exhalation resistance 160 l/min (mbar)		<3 mbar	<3 mbar	<3 mbar

PPE's controlled by:

BSI Group The Netherlands B.V. (2797), Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, Netherlands

Our filtering half masks are in conformity with the PPE Regulation(EU) 2016/425.

EU Declaration of Conformity is available to download from our website under:

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**EN149:2001+A1:2009**

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**FFP1 NR D 4910**  
**FFP2 NR D 4920/4921**  
**FFP3 NR D 4931**

SUITABLE FOR USE AGAINST SOLID,WATER AND LIQUID BASED AEROSOLS.

**INSTRUCTIONS FOR USE:** Single-use filtering face-piece respirators

**IMPORTANT:** Single-use respirators are classified according to one of three categories-FFP1 FFP2, FFP3. Before reading the following information, check the category to which the respirator belongs-this is indicated on the packaging and on the respirator.

If the mask is marked D, it has passed the dolomite clogging test and suitable for use in very dusty environment.

It is the user's responsibility to ensure that the respirator provides the necessary level of protection for the type and concentration of the contaminant(s) in the area where the respirator is intended for use.

#### DIRECTIONS FOR USE:

- Fit the respirator and check leak-tightness BEFORE entering the contaminated area.
- Wear the respirator for the whole duration of exposure to contaminants.
- Use the respirator in accordance with applicable health and safety regulations.
- Discard the respirator and replace with a new one if: \*the respirator is removed while in a contaminated area \*excessive clogging of the respirator causes breathing difficulty or discomfort \*the respirator becomes damaged \*(for respirators protecting against vapours)-the smell of vapours present becomes detectable.
- Leave the contaminated area if dizziness, irritation or other distress occurs.
- Only for single use. No maintenance necessary. Do not store and/or re-use after single use.
- Discard the respirator after single use (one shift, maximum).
- Keep un-used respirators in their closed box and store in a dry non-contaminated area between -20 and +40 deg. C. at relative humidity under 75%

NOTIFIED BODY: BSI Group The Netherlands B.V. (2797), Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, Netherlands

#### FITTING INSTRUCTIONS:

1. Cup the respirator in the hand, allowing the head-bands to hang freely.
2. Hold the respirator under chin with the nose-piece facing outwards.
3. Place the lower head-band around the neck below the ears. Holding the respirator against the face with one hand, place the top head-band above the ears, around the crown of the head.
4. Respirators with adjustable head-bands only: While maintaining the position of the respirator, adjust tension by pulling tab of each head-band. (Tension may be decreased by pushing out on the back of the buckle).
5. Mould the nose area to the shape of the face, running the fingertips of both hands from the top of the nose-piece down both sides while pressing inward.

**CHECK FOR LEAK-TIGHTNESS (FACIAL FIT) AS FOLLOWS:** \*place both hands over the respirator and exhale sharply \*positive pressure should be inside the respirator \*if you feel air escaping around the edges, re-adjust the respirator by lightening the nose-piece and/or increasing tension of the head-bands \*repeat the procedure until satisfactory face-fit has been achieved.

**IMPORTANT:** In the case of respirators with non-adjustable head-bands-if you feel that the head-bands are too loose, you can increase tension by tying a small knot in the head-band.

#### LIMITATIONS TO USE:

- Do not use these respirators or enter or stay in an area where:
- Oxygen concentration is less than 17% (UK-19%)
- Contaminants or their concentrations are unknown or immediately dangerous to life or health
- Particulate concentrations exceed levels fixed by applicable health and safety regulations or protection factor x NPF-whichever is lower. (Nominal protection factors: 4.5 for FFP1, 12.5 for FFP2 and, 50 for FFP3 )
- Gases and/or vapours are present-except for respirators designated as protecting against vapours, in which case vapour concentrations should not exceed NPF. Respirators protecting against organic vapours should not be used as protection against acid vapours and vice-versa.

#### WARNING:

Do not use for fire-fighting. These respirators do not supply oxygen. Do not use in oxygen-deficient atmospheres-e.g. tanks or other poorly-ventilated areas (see Limitations to Use). Do not use in explosive atmospheres. For use only by trained and qualified personnel. Requirements for leak-tightness are unlikely to be achieved if the respirator is worn against a beard or facial stubble.

When transporting this product, retain in the original packaging and keep away from mechanical and chemical hazards.